DEUTSCH

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






ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE:

The intended user should be a licensed medical professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

INDICATIONS: The DonJoy Telescoping Post-Op Brace is designed to aid in immobilization of the leg and provide protected range of motion following ACL, PCL, LCL and MCL surgeries and meniscal repairs. Providing immobilization or controlled movement of the limb or body segment. Providing mild protection of a limb or body segment for acute and prophylactic care.

CONTRAINDICATIONS: None

APPLICATION INSTRUCTIONS:

- 1) Unfasten the buckle on each strap of the brace and detach hinge bars from foam liners. Attach the hook end of the foam liners to the lateral hinge bar with a 4" (9 cm) tail of foam liner on the front.
- 2) Adjust bar length by pushing in the gray button and sliding the extension bar to desired length. The button must fully extend up after adjustment to insure a secure position and both bars must be equal in length. The post-op bars may be extended the full length of the leg to the ankle.
- 3) Position the lateral side (outside) bar assembly aligning the hinge with the patella (kneecap) and midline of the leg. The lateral malleolus and greater trochanter at the hip serve as good reference points for midline alignment.
- 4) Wrap foam calf and thigh pads snugly around leg and attached anteriorly.
- 5) Position and press into place the medial (inside) bar. Make sure medial and lateral hinge height is equal.
- 6) Adjust straps lengths to the leg, and then snap-in buckles starting with the two straps closest to the knee. Cinch straps to secure in place. Fasten remaining buckles and straps.
- 7) For easy one-piece removal, unfasten the buckles and open the foam pads. The brace can then be removed and reapplied as a single unit.

INSTRUCTIONS FOR BAR ADJUSTMENT

The length of the thigh and calf bars may be adjusted to fit the length of the leg. Adjust bar length by pushing in the gray button and sliding the extension bar to desired length. The button must fully extend up after adjustment to insure a secure position and both bars must be equal in length. Bars may need to be readjusted for malleolus comfort. The extension bar may be removed by carefully sliding the bar completely out.

IROM HINGE: A Press down, slide out, and lift up the hinge security latch. Rotate the cover until the desired flexion or extension pin is exposed. Remove the pin and place it in the desired degree setting. Rotate the cover so that the security latch is past the pin and replace the latch.

ELS HINGE: B To engage the hinge, the leg must be fully extended. Push the button on the back of the hinge inward and then forward. To unlock the hinge, push the button on the front of the hinge inward and then back.

TROM HINGE: C To adjust flexion and extension, push in button then adjust tab to desired position. The hinge bar may need to be moved to allow adjustment. To lock hinge, adjust to desired position, push in button and insert tab into lock position. NOTE: The button must fully extend up after adjusting the flexion and extension angle or inserting the lock to insure a secured position.

USE AND CARE

- The hinge bars may be bent to add varus or valgus contouring. Bend each side bar an equal amount above and below the hinge.
- Hand wash foam liners in cold water with mild detergent. AIR DRY ONLY, DO NOT HEAT DRY.
- All foam may be trimmed.

WARNINGS AND PRECAUTIONS:

- Loss of circulation, patient discomfort and patient re-injury are potential effects caused by device failure.
 - ROM settings should not be changed without supervision of a medical professional.
 - This product must be prescribed and fitted by a healthcare professional.
 - The frequency, duration of use and directions for use should be determined by your healthcare professional.
 - If you develop an allergic reaction and/or experience itchy, red skin after coming into contact with any part of this device, please stop using it and contact your healthcare professional immediately.
 - If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.
 - The support should be snug but not impair circulation.
 - Do not use over open wounds.
 - Do not use the device if it is damaged and/or the packaging has been opened.
- NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

INTENDED FOR USE ON A SINGLE PATIENT.

RX PRODUCT

NOT MADE WITH NATURAL RUBBER LATEX.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THIS DEVICE IS ONLY ONE ELEMENT IN THE OVERALL TREATMENT PROGRAM ADMINISTERED BY A MEDICAL PROFESSIONAL. THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PROFESSIONAL.

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©2021 DJO, LLC
MDSS GmbH