



## ENGLISH

**BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.**

### INTENDED USER PROFILE:

The intended user should be a licensed medical professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

### INTENDED USE/INDICATIONS:

The DonJoy Telescoping TROM Advance Knee Brace is designed to aid in immobilization of the knee and provide protected range of motion following ACL, PCL, LCL and MCL surgeries and meniscal repairs. Providing immobilization or controlled movement of a limb or body segment. Providing mild protection to a limb or body segment for acute and prophylactic care.

### CONTRAINDICATIONS:

None

### WARNINGS AND PRECAUTIONS:

- Loss of circulation, patient discomfort and patient re-injury are potential effects caused by device failure.
- ROM settings should not be changed without supervision of a medical professional.
- This product must be prescribed and fitted by a healthcare professional.
- The frequency, duration of use and directions for use should be determined by your healthcare professional.
- If you develop an allergic reaction and/or experience itchy, red skin after coming into contact with any part of this device, please stop using it and contact your healthcare professional immediately.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.
- The support should be snug but not impair circulation.
- Do not use over open wounds.
- Do not use the device if it is damaged and/or the packaging has been opened.

NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

### APPLICATION INFORMATION:

- Unfasten the buckle on each strap of the brace and remove foam wraps from brace. Wrap foam calf and thigh pads snugly around leg and attach anteriorly (front of leg). The foam can be trimmed for proper fit.
- Adjust bar length by pushing in the gray button and sliding the extension bar to desired length. The button must fully extend up after adjustment to ensure a secure fit. Both bars must be equal in length. The post-top bars may be extended the full length of the leg to the ankle.
- Position and press into place the lateral side (outside) bar assembly, aligning the hinge with the patella (kneecap) and midline of the leg. The lateral malleolus and greater trochanter at the hip serve as good reference points for midline alignment. Repeat for the medial (inside) bar. Make sure medial and lateral hinge height is equal.
- To adjust the position of the straps on the leg, push the gray button and slide the Tele-Fit strap paddles to the desired location on the leg.
- Fasten the two straps closest to the knee, then proceed with the remaining two straps. Loosen the straps by pulling the straps away from the frame, and pull straps tight to remove slack behind the leg while maintaining the position of the hinge and bars on the leg. The straps can be trimmed for proper fit and strap-end can be repositioned to new edge of strap.
- To adjust flexion and extension, slide the buttons found on the hinge casing to the desired flexion/extension settings as determined by your prescribing healthcare professional. You may need to move the hinge bar to the desired extension setting on the outer rim of the hinge in order to be able to move the extension button.
- To lock the hinge, move the hinge bar to the desired setting on the outer rim of the hinge. Slide the hinge lock button, located on the thigh bar, to the locked position such that the arrow moves behind the hinge center.
- For easy one-piece removal, unfasten the buckles and open the foam pads. The brace can then be removed and reapplied as a single unit.
- The hinge may be bent to add or subtract valgus/corriging. Bend each side bar an equal amount above and below the hinge.
- For optimal ease of application post-operatively, pre-fit the brace prior to surgery if possible.

**CLEANING INSTRUCTIONS:** Hand wash foam liners in cold water with mild detergent. Air dry only, do not heat dry. Regular cleaning of the brace is recommended. All foam may be trimmed.

**WARRANTY:** DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

### INTENDED FOR USE ON A SINGLE PATIENT.

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**NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN AN EFFORT TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.**

**CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PROFESSIONAL.**

**NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN AN EFFORT TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.**

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