

enovis™

DONJOY®
SUPPORT EVEREST® II

Begin your ascent



SUPPORT EVEREST™ II

Support Everest™ II from DonJoy® offers an enhanced option for support and protected range of motion (ROM) control of the knee.

INTENDED USE/INDICATIONS

The DonJoy Support Everest II is designed to support the knee and provide protected ROM following light to moderate ACL/PCL injuries, moderate to severe lateral ligament injuries (MCL, LCL), and meniscal injuries.

FEATURES AND BENEFITS

• Improved polycentric bilateral hinges

- For increased ROM control with tool-free flexion/extension stops at 0°, 10°, 20°, 30°, 60° and 90°. Full ROM also available.
- Bilateral stays help provide medial / lateral stabilization

• Improved comfort

- Breathable elasticated mesh fabric promotes airflow to help provide comfort
- Lightweight, slimline material also helps reduce weight and profile of brace

• Improved fit

- Contoured anatomical fit around the patella
- Available in 6 sizes (XS to XXL), fits left or right leg
- 2 different brace heights: XS/S/M = 39.5 cm, L/XL/XXL = 44.5 cm

• Removable posterior stay

- Rigid enough for strict immobilization, malleable enough to be formed to patient's needs
- Can be removed in rehab phase, to assist un-locking of knee joint

• Improved buckles

- Easy to fasten and adjust

• Wraparound structure

- Easy to apply and adapts to patient anatomy

• 4 independent straps

- To adjust level of compression and support

• Anti-slip strap


- Internal band positioned under the knee helps maintain correct position



ORDERING INFORMATION

PART NUMBER	THIGH CIRCUMFERENCE		SIZE
	CM	Inches	
82-0046-1	33 – 37 cm	13 – 14½"	XS
82-0046-2	38 – 43 cm	15 – 17"	S
82-0046-3	44 – 49 cm	17¼ – 19¼"	M
82-0046-4	50 – 55 cm	19½ – 21½"	L
82-0046-5	56 – 60 cm	22 – 23½"	XL
82-0046-6	61 – 64 cm	24 – 25¼"	XXL

Measurement taken 15 cm above mid-patella



enovis™

Creating Better Together™

enovis-medtech.eu

2000 Cathedral Square, Cathedral Hill, Guildford, Surrey, GU2 7YL, UK

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.