



Exos FORM™ Spine Bracing System

Reshaping Form, Fit and Function

exos
form

Shaping a Revolution

Innovation fueled by passion, guided by the functional beauty of the human form. Refined simplicity with smooth, beautiful lines. An optimized fit tuned to each body. Compromise nothing. Challenge the status quo. Reinvent, redefine, reimagine. Above all, form honors function. These were the inspiration for a visionary shift in spine bracing. Exos™ is unrelenting in its quest to find new ways to shape care, comfort and healing in the human body, and the Exos™ Spine Brace System is an opus in engineering, craftsmanship and detailed precision, built to deliver an unparalleled progression of modular support.



Exos™ Thermoformable Panels

A quick and easy way to ensure a personalized fit.

Simple, Light Weight, Low Profile Design

Designed with patient comfort in mind without sacrificing functionality.



Independent Superior & Inferior Compression

Superior spinal decompression delivered through Boa® Technology's innovative and elegant closure system.

Adjustable (Angle) Belt Wings

The ability to angle the belt to accommodate different waist/torso contours to ensure an optimal fit.


Modular Step-Up/Step-Down Design

The ability to adjust and configure your brace as you progress through rehabilitation and recovery.

Semi-Universal Design

Minimizes inventory management requirements.





Exos™ Spine Bracing System

The Exos™ Spine Brace transforms spine bracing. Designed with Exos™ thermoformable technology and Boa® Technology's advanced closure system, the Exos™ Spine Brace combines superior decompression and support, with the ability to quickly and easily customize fit. The end result is a spine brace that will change the way you deliver patient comfort and satisfaction.



Exos FORM™ 627

The Exos FORM™ 627 (with anterior panel) provides relief from a wide range of indications from mild to acute lower back pain, as well as post-operative support. The Exos FORM™ 627 provides superior sagittal control and support from L-1 to L-5.

PART NO.	DESCRIPTION
205627-40	EXOS FORM™ 627, S/M
205627-60	EXOS FORM™ 627, L/XL
PDAC Assigned Code: Lo627	



Exos FORM™ 637

The Exos FORM™ 637 provides a higher degree of spinal support and relief from a wide range of indications from acute lower back pain to post-operative support. The Exos FORM™ 637 provides superior sagittal and coronal control and support from T9-S1.

PART NO.	DESCRIPTION
205637-40	EXOS FORM™ 637, S/M
205637-60	EXOS FORM™ 637, L/XL
PDAC Assigned Code: Lo637	



Exos FORM™ 626

The Exos FORM™ 626 that provides relief from mild lower back pain. The Exos FORM™ 626 provides superior sagittal control and support from L-1 to L-5.

PART NO.	DESCRIPTION
205626-40	EXOS FORM™ 626, S/M
205626-60	EXOS FORM™ 626, L/XL
PDAC Assigned Code: Lo626	



Exos FORM™ 631

The Exos FORM™ 631 provides a higher degree of spinal support and relief from a wide range of indications from acute lower back pain to post-operative support. The Exos FORM™ 631 provides superior sagittal control and support from T9-S1.

PART NO.	DESCRIPTION
205631-40	EXOS FORM™ 631, S/M
205631-60	EXOS FORM™ 631, L/XL
PDAC Assigned Code: Lo631	



Exos FORM™ 621

The Exos FORM™ 621 provides superior stabilization of the sacroiliac (SI) joint through direct compression. Silicon grip strips prevent the brace from migrating.


PART NO.	DESCRIPTION
205621-40	EXOS FORM™ 621, S/M
205621-60	EXOS FORM™ 621, L/XL
Recommended PDAC Code: Lo621	

COMMON INDICATIONS

INDICATIONS

	Sacrum	L1 - L5	L1 - L5	T9 - S1	T9 - S1
	Exos FORM™ 621	Exos FORM™ 626	Exos FORM™ 627	Exos FORM™ 631	Exos FORM™ 637
Sacroiliac Pain	●				
Acute and Chronic Lower Back Pain		●	●	●	●
Lower Back Sprains / Strains		●	●	●	●
Lumbar Disc Displacement		●	●	●	●
Osteoporosis		●	●	●	●
Disc Herniation and Degeneration		●	●	●	●
Post-Operative Laminectomy		●	●	●	●
Post-Operative Disectomy		●	●	●	●
Spondylolisthesis		●	●	●	●
Post-Operative Fusion				●	●
Spinal Stenosis				●	●
Stable / Non-Displaced Fractures				●	●
Compression Fractures					
Scoliosis					

RANGE OF SPINAL SUPPORT





ORDERING INFORMATION

SIZING GUIDE WAIST CIRCUMFERENCE

DESCRIPTION	Exos FORM™ 637	Exos FORM™ 631	Exos FORM™ 627	Exos FORM™ 626	Exos FORM™ 621
S/M	25" - 49" (63cm - 125cm)	25" - 49" (63cm - 125cm)	25" - 49" (63cm - 125cm)	25" - 49" (63cm - 125cm)	less than 35" (89 cm)
L/XL	49" - 60" (125cm - 152cm)	49" - 60" (125cm - 152cm)	49" - 60" (125cm - 152cm)	49" - 60" (125cm - 152cm)	36" - 55" (91 cm - 140 cm)

EXOS FORM™ SPINE BRACING SYSTEM

PART NO.	DESCRIPTION	COLOR	SIZE
205637-40	EXOS FORM™ 637, S/M	BLACK	S/M
205637-60	EXOS FORM™ 637, L/XL	BLACK	L/XL
205631-40	EXOS FORM™ 631, S/M	BLACK	S/M
205631-60	EXOS FORM™ 631, L/XL	BLACK	L/XL
205627-40	EXOS FORM™ 627, S/M	BLACK	S/M
205627-60	EXOS FORM™ 627, L/XL	BLACK	L/XL
205626-40	EXOS FORM™ 626, S/M	BLACK	S/M
205626-60	EXOS FORM™ 626, L/XL	BLACK	L/XL
205621-40	EXOS FORM™ 621, S/M	BLACK	S/M
205621-60	EXOS FORM™ 621, L/XL	BLACK	L/XL

EXOS FORM™ SPINE BRACING ACCESSORIES

PART NO.	DESCRIPTION	COLOR	SIZE
805-631-40	POSTERIOR T9 PANEL, S/M	BLACK	S/M
805-631-60	POSTERIOR T9 PANEL, L/XL	BLACK	L/XL
805-637-40	LATERAL EXTENSIONS, S/M	BLACK	S/M
805-637-60	LATERAL EXTENSIONS, L/XL	BLACK	L/XL
805-AP-40	ANTERIOR PANEL, S/M	BLACK	S/M
805-AP-60	ANTERIOR PANEL, L/XL	BLACK	L/XL
805-1-WINGS-40	BELT WING SET, S/M	BLACK	S/M
805-1-WINGS-60	BELT WING SET, L/XL	BLACK	L/XL
805-1-BWING-40	BRIDGE WING, S/M	BLACK	S/M
805-1-BWING-60	BRIDGE WING, L/XL	BLACK	L/XL
805-1-INSERT-40	LORDOTIC INSERTS, S/M	BLACK	S/M
805-1-INSERT-60	LORDOTIC INSERTS, L/XL	BLACK	L/XL
805-JITT	SI BELT PAD SET	BLACK	UNIVERSAL

Embrace the Reform



DJO, LLC | A DJO Global Company

T 800.336.6569 **D** 760.727.1280 **F** 800.936.6569
1430 Decision Street | Vista, CA 92081-8553 | U.S.A.

DJOglobal.com/exos

Together in Motion™



U.S.A. and International Patents Pending

©2013 DJO, LLC

00-1262 Rev A