

enovis™

DONJOY®

DONJOY IMMOSTRAP 2.0™

Firmer, targeted
lumbar support



DONJOY®

DONJOY IMMOSTRAP 2.0™

Firmer, targeted lumbar support

Active patients with low back pain or a herniated disc may benefit from using DonJoy Immostrap 2.0™. Thanks to its independent upper and lower bilateral straps and flexible dorsal stays, this highly adjustable lumbar support can help deliver pain relief and firmer targeted postural control.

INTENDED USE/INDICATIONS

The DonJoy Immostrap 2.0 is designed to provide comfort and support of the lumbar spine. It may be suitable for the symptomatic relief of acute and chronic low back pain. Provides compression without significant restriction of mobility, can be adjusted as and when required by the adjustable strap.

- Acute lower back pain
- Chronic low back pain
- Herniated disc

FEATURES AND BENEFITS

- **Now available in 2 heights 26 cm and NEW 21 cm**, allowing enhanced fit for male and female patients
- **Dual hand loops** for easy application
- **4 flexible dorsal stays** provide support to the lower back
- **Adjustable double bilateral straps** for firmer postural control



HEIGHT 21 cm		HEIGHT 26 cm		WAIST CIRCUMFERENCE	SIZE
PART NUMBER	GTIN	PART NUMBER	GTIN		
82-0074-1	0190446825234	82-0075-1	0190446826514	56 - 68 cm	XS
82-0074-2	0190446826460	82-0075-2	0190446826521	68 - 82 cm	S
82-0074-3	0190446826477	82-0075-3	0190446826538	82 - 98 cm	M
82-0074-4	0190446826484	82-0075-4	0190446826545	98 - 116 cm	L
82-0074-5	0190446826491	82-0075-5	0190446826552	116 - 136 cm	XL
82-0074-6	0190446826507	82-0075-6	0190446826569	136 - 156 cm	XXL

