

# AirSport+™

*Simple. Stable. Smarter.*

---

The **Aircast® AirSport+** is a clinically proven stirrup design providing ankle stabilization for moderate to severe ankle sprains and chronic instabilities. The anatomically correct design and durable, breathable materials offer a stable, comfortable fit. Along with the Aircast technology, simple step in design and intuitive, unified closure system, the AirSport+ provides a functional management option to address ankle injuries and rehabilitation.

- ◆ **Low profile, customizable fit** with BOA® Lacing System.
- ◆ Anterior TaloFibular (ATF) **cross strap provides compression** and stabilization to the ankle joint, helping to minimize inversion.
- ◆ **Anatomically designed** medial/lateral semi-rigid shells.
- ◆ Posterior high ankle strap – provides **compression to the tibia/fibula**.
- ◆ **Simple, easy application** with posterior step in design.



# AirSport+

*Simple. Stable. Smarter.*



			SHOE SIZE	
RIGHT	LEFT	SIZE	MALE	FEMALE
02RXSR	02RXSL	XS	< 3	< 5
02RSR	02RSL	S	3 – 5	5 – 6 ½
02RMR	02RML	M	5 – 8 ½	6 ½ – 10
02RLR	02RLL	L	8 ½ – 12	10 – 13 ½
02RXLR	02RXLL	XL	>12	>13 ½

**For more information on the new AirSport +,  
please contact your local DJO Territory Manager**



**DJO, LLC | A DJO Global Company**

**T** 800.336.6569 **D** 760.727.1280 **F** 800.936.6569

1430 Decision Street | Vista, CA 92081-8553 | U.S.A.

**[DJOglobal.com](http://DJOglobal.com)**