DONJOY®

MyBabyStrap®

Enjoy every day of your pregnancy





MyBabyStrap®

Lumbar maternity belt to aid in relief and support during pregnancy and after birth

Indications: Pain and discomfort associated with pregnancy

- \cdot Lower back pain
- $\cdot \, \text{Pelvic pain}$
- \cdot Sacroiliac pain
- \cdot Postural instabilities
- Discomfort due to stretching and expansion of the abdomen
- Back support after birth



Modular design:

adaptable to all stages of pregnancy and also after birth. Easy to adjust

Adjustable back pad:

to help provide additional support and compression in the back area



Four posterior rigid stays: for extra support in the

lumbar area

Bilateral straps:

to help provide adjustable relief. Tightness can be adjusted according to individual needs and level of pain



Anatomical shape: follows the natural contours of a pregnant body for an improved fit, with no pressure on the abdomen





Stylish maternity belt: featuring feminine details (silver threads and scalloped edge design)



Thin, soft and breathable fabric: for a comfortable fit



Ergonomic front panel: for comfort when sitting or standing



Easy-to-use with simple, hook-loop closure: patients can quickly adjust the belt to the desired level of comfort



Removable elastic strap: for lower abdominal support. Easy application



DONJOY®



MyBabyStrap®

Part Number	Waist Circumference
82-0276	80 – 150 cm

Ventral height: 8 cm Dorsal height: 22 cm Machine washable



DJO Global | 1a Guildford Business Park | Guildford | Surrey | GU2 8XG | UNITED KINGDOM

www.DJOglobal.eu